


Ready, set, GO!

Transform your EMS studio with our cutting-edge solutions. Beyond just weight, our products emphasize the positive effects on body composition through EMS training. When clients witness tangible progress and improvement, they're more inclined to increase session bookings and refer others.

Improve EMS training experience with the **TANITA MC-780** segmental body composition monitor. This device offers a **medical-grade analysis** of clients' health and fitness, providing insights beyond mere weight measurement.

With **20+ comprehensive measurements** and tracking progress, you can diagnose areas for improvement and guide clients towards achieving their goals effectively.



JAPANESE TECHNOLOGY | **5 YEAR GUARANTEE** | **WORLD NO. 1**
FOR NEW HEALTH INSIGHTS

Key Features

- **Full segmental body composition measurements** performed in under 20 seconds using the latest multi frequency BIA technology
- **Extra large LED dual display** show analysis in an easy-to-read format.
- **Data can be sent wirelessly** to a PC or Windows tablet using a Bluetooth adaptor.
- **Max weight capacity** 270kg with 100g accuracy
- **Results automatically stored** on the SD Card, sent to a PC or printed directly with a PictBridge printer as a complete consultation sheet for further discussion.
- **Interactive stand-alone unit** allows clients to quickly take a measurement without specialist assistance
- **Accuracy grade:** NAWI Class III, MDD Class II-a

Measurements*

- Weight
- Body Mass Index (BMI)
- Body Fat (kg / %)
- Visceral Fat Indicator
- Fat Mass
- Fat Free Mass
- Muscle Mass
- Protein (kg)
- Total Body Water (kg / %)
- Extra-Cellular Water (kg)
- Intra-Cellular Water (kg)
- ECW/TBW Ratio
- Basal Metabolic Rate
- Basal Metabolic Rate Indicator
- Bone Mineral Mass Indicator
- Metabolic Age
- Physique Rating
- Phase Angle
- Sarcopenia Index

Segmental Measurements

- Body Fat (kg / %)
- Fat Distribution Rating
- Muscle Mass (kg)
- Muscle Mass Rating
- Muscle Mass Balance
- Leg Muscle Score
- Reactance/Resistance
- Phase Angle

* Display of measurements may vary: Printer, TANITA PRO software or TANITA PRO app may be required.



Also available in a portable format

MC-780MA S

Instead, the easy-to-read display can be placed on your desk for example, while the scale can remain on the floor.

